



Invacare eForms offer new functionality!

Invacare is pleased to offer you an upgraded solution to your ordering process. Our enhanced order forms allow you to fill out a form electronically, print and fax the form, save and email* it to Customer Service, or maintain the business practices that work for you today. The format has been revised to reveal a cleaner look with electronic selection and input functions.

*Saving and emailing the eForms is available with Acrobat Professional software

If you are using Adobe Reader

- **View** and download all forms
- **Save** blank versions of the forms to your desktop
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service

Please note: The 'Email' button at the end of the form will not work with Adobe Reader. The 'Email' button will only work while using Adobe Acrobat Professional.

If you are using Adobe Acrobat Professional

- **View** and download all forms
- **Save** blank AND completed versions of the forms to maintain electronic records
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service
- **Email** completed forms to Customer Service by using the 'Email' button at the end of the form.

For information on downloading Adobe Reader for free, or purchasing Adobe Acrobat Professional, please visit www.adobe.com.

Invacare® 3G Arrow® Power Wheelchair with TrueTrack

Price List and Order Form

www.invacare.com

Customer Service: 800-333-6900

Fax: 800-678-4682

To order literature: 800-828-6282



TRANSACTION TYPE:

Order Quote

Date _____ Account # _____ P.O. # _____

Company Name _____ Phone # _____

Address _____ City _____ State _____ Zip Code _____ Country _____

Contact Name _____ Comments _____

Special Note

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.

CHAIR TYPE

- 3GAR Arrow Base and Seat with MK6i™ \$9,775.00 ^{K0868}
Arrow comes standard with 300 lbs. or 400 lbs. user weight capacity, 7 MPH top speed, TrueTrack driving technology with GB™ motors, MK6i CMPJ™+ joystick, MK6TT controller, Base only does not come standard with Seat Positioning Strap.
- 3GARBASE Arrow Base Only with MK6i \$9,020.00
Arrow comes standard with 300 lbs. or 400 lbs. user weight capacity, 7 MPH top speed, TrueTrack driving technology with GB™ motors, MK6i CMPJ™+ joystick, MK6TT controller, Base only does not come standard with Seat Positioning Strap.
- 3GSEAT Seat Only \$755.00
When ordering a replacement seat and changing depth to 20" or deeper, must also order new arms.

TRANSPORT TIE-DOWN Clear

- TRRO Transport Ready Option \$300.00
- TRBKTS Wheelchair Transport Brackets \$250.00

Group 4 - K0868 Basic Package	
3GAR Arrow with MK6i	
300 lbs. Weight Capacity	
1431-3 14" x 3" Tires with Foam Filled Inserts	No Charge
TrueTrack Gearless Brushless GB™ Motors	Standard
1228-3 8" x 2 1/4" Semi Pneumatic Casters	Standard
28F Desk Length Fixed Height Arm	Standard
MK6TT TrueTrack Controller	Standard
70NHD 70° Swingaway Footrests	Standard
8 Amp Off Board Battery Charger	Standard
1515 Push Button Style Seat Positioning Strap	Standard

Group 4 - K0870 Basic Package	
3GAR Arrow with MK6i	
400 lbs. Weight Capacity	
1431-3 14" x 3" Tires with Foam Filled Inserts	No Charge
TrueTrack Gearless Brushless GB™ Motors	Standard
1228-3 8" x 2 1/4" Semi Pneumatic Casters	Standard
28F Desk Length Fixed Height Arm	Standard
MK6TT TrueTrack Controller	Standard
70NHD 70° Swingaway Footrests	Standard
8 Amp Off Board Battery Charger	Standard
1515 Push Button Style Seat Positioning Strap	Standard

Transport Tie Down*

TRRO includes four factory-installed transport brackets and a wheelchair anchored pelvic belt. TRRO has been crash-tested in accordance with ANSI/RESNA WC Vol 1. Section 19 Frontal Impact Test requirements for wheelchairs with a 168 lbs. crash dummy, which corresponds to a person with a weight of 114 lbs. to 209 lbs.

TRBKTS includes four factory-installed wheelchair transport brackets. TRBKTS has not been crash-tested in accordance with WC19. Use these transport brackets only to secure an unoccupied wheelchair during transport.

As of this date, the Department of Transportation has not approved any tie-down systems for transportation of a user while in a wheelchair, in a moving vehicle of any type. It is Invacare's position that users of wheelchairs should be transferred into appropriate seating in vehicles for transportation and use be made of the restraints made available by the auto industry. Invacare cannot and does not recommend any wheelchair transportation systems.

* Transport Tie-Down option is not retrofittable to existing models and is not field serviceable.

TRRO Features for 3G Arrow

- | Required Features for TRRO | Optional Features for TRRO |
|--|---|
| • Storm Series Arrow Base & Seat | • Fender |
| • 300 lb or 400 lb user weight limit | • Flip-Up Composite Footplates, Aluminum Footplates-Deluxe or Adjustable Angle Flip-Up Footplates |
| • Short or Long Base Frame | • Calf Strap |
| • Adjustable Rehab Seat with Adjustable Seat Tilt & Back Angle | • Impact Guards |
| • Low Seat-to-Floor Height | • All PSR/PSF Joysticks |
| • Sectional or Straight Rear Uprights | • Left Hand Mounted Joystick |
| • 18"- 24" ASBA Back Angle | • "T" Handle Flexible Joystick Ext. |
| • 8" x 2.25" Semi Pneumatic Casters | • Straight Handle Flexible Joystick Ext. |
| • Reinforced Back Uph. | • MK6i Remote Programmer |
| • 14" x 3" Tires w/Foam Filled Inserts | • Group 24 Batteries |
| • Rear Wheel Standard Location | |
| • Black Nylon or Black Vinyl Uph. | |
| • Center Pivot Style Bracket | |
| • 60° Swingaway, 70° Swingaway or 70° Tapered Footrests | |

Please note that all TRRO chairs are set at 10° seat angle

**NOTE: All specifications and dimensions are approximate.
 For the most current pricing information, go to www.invacare.com/homecare.**

Invacare® 3G Arrow® Power Wheelchair w/TrueTrack

FRAME SELECTION

USER WEIGHT LIMITS

- U300 User Weight < 300 lbs. No Charge K0868
- U400 User Weight < 400 lbs. No Charge K0870

ADJUSTABLE SEAT WIDTH RANGE

- Clear
- 1216W 12" to 16" Seat Width Range No Charge
- 1620W 16" to 20" Seat Width Range No Charge
- 2024W 20" to 24" Seat Width Range No Charge

FRAME TYPE

- Clear
- SFS Short Frame No Charge
SFS Not Available on seat depths of 1922D (19" to 22" seat depth range).
- LFS Long Frame No Charge

FRAME WIDTH

- Clear
- For Frame Widths of 12"-15":**
Jr. Seat widths/depths require a straight rear upright and a user weight limit of 150 lbs. Manual Reclining back not available. Additional charges apply to widths/depths less than 16". Jr. seat widths have Jr. arms with an adjustment range of 6"-10" inches in height.
- W12 12" Wide \$110.00
- W13 13" Wide \$110.00
- W14 14" Wide \$110.00
- W15 15" Wide No Charge
- W16 16" Wide No Charge
- W17 17" Wide No Charge
- W18 18" Wide No Charge
- W19 19" Wide No Charge
- W20 20" Wide No Charge
- W21 21" Wide \$329.00
- W22 22" Wide \$329.00
- W23 23" Wide \$329.00
The 23" and 24" seat widths are not available with low seat-to-floor. Must pick Heavy-Duty Package to order 23" and 24" seat widths.
- W24 24" Wide \$469.00
The 23" and 24" seat widths are not available with low seat-to-floor. Must pick Heavy-Duty Package to order 23" and 24" seat widths.

SEAT SELECTION

SEAT-TO-FLOOR HEIGHT

Listed values are for 0 degree Seat Angle and w/o TRBKTS. Please see below for additional information. (Values with an "F" represent factory settings, which are the lowest possible STF)

- LOW Seat-To-Floor - Low 17.25" No Charge
- MED Seat-To-Floor - Med 19.50" No Charge
- TALL Seat-To-Floor - Tall 20.75" No Charge

SEAT SELECTION

SEAT SELECTION

- ADJASBA Adj Rehab Seat with Adj Seat Tilt and Back Angle Standard
- ADJRECL Manual Recline for Adjustable ASBA \$744.00 E1226
Not available with elevating seat.
- VAN Van Style Captain's Seat \$100.00
Full length adjustable arms are standard on Van Style Seat.
- CURV Cruved Back No Charge E1226
- CTB Contoura Seating..... \$625.00 E2620 E2621
 - OMCU Omit Seat Cushion Less \$20.00
- CTBMSH Contoura seating with mesh upholstery \$635.00 E2620 E2621
 - OMCU Omit Seat Cushion Less \$20.00

VAN SEAT STYLE

- Clear
- Used only on Van Style seat selection.**
- SEMIR Semi-Recline No Charge
- FULLR Full-Recline \$469.00 E1226

VAN SEAT WIDTH

- Clear
- Manual Reclining Back not available.**
Van Seat Depth set at the following:
18" wide set @ 16" deep
20" wide set @ 18" deep
22" wide set @ 20" deep
- VW18 18" Wide No Charge
- VW20 20" Wide \$329.00
- VW22 22" Wide \$329.00

ADJUSTABLE SEAT DEPTH RANGE

- Clear
- 1218D 12" to 18" Seat Depth Range No Charge
Must select 12" - 18" depth ranged (1218D) when ordering the 12" - 16" seat width range (1216W). 12" and 13" widths are not available with 20" - 24" back heights.
- 1619D 16" to 19" Seat Depth Range No Charge
- 1922D 19" to 22" Seat Depth Range No Charge

SEAT DEPTH SETTING

- Clear
- Contoura Seating available in depths 16" - 22".**
For Seat Depths of 12"-15":
Jr. Seat widths/depths require a straight rear upright and a user weight limit of 150 lbs. Manual Reclining back not available. Additional charges apply to widths/depths less than 16". Jr. seat widths have Jr. arms with an adjustment range of 6" - 10" inches in height.
- D12 Seat Depth 12" (Short Base) \$110.00
- D13 Seat Depth 13" (Short Base) \$110.00
- D14 Seat Depth 14" (Short Base) \$110.00
- D15 Seat Depth 15" (Short Base) No Charge
- D16 Seat Depth 16" (Short Base) No Charge
- D17 Seat Depth 17" (Short Base) No Charge
- D18 Seat Depth 18" (Short Base) No Charge
- D19 Seat Depth 19" (Long Base) No Charge
- D20 Seat Depth 20" (Long Base) No Charge
- D21 Seat Depth 21" (Long Base) \$529.00
- D22 Seat Depth 22" (Long Base) \$529.00

		Short Frame (12-16" or 16-19" Seat Depth)			Long Frame (19-22" Seat Depth)			
		Low	Medium	Tall	Low	Medium	Tall	
Non - Powered Seating	0 Degree Seat Angle	w/ TRBKTS	17.5	19.75F	21F	17.75	20F	21.25F
		w/o TRBKTS	17.25F	19.5F	20.75F	17.5F	19.75F	21F
	5 Degree Seat Angle	w/ TRBKTS	18.25F-19	20.5-21.25	21.75-22.5	18.75F-19.5	21-21.75	22.25-23
		w/o TRBKTS	18-18.75	20.25-21	21.5-22.25	18.5-19.25	20.75-21.5	22-22.75
		w/ TRRO	19.50F			20F		

For the most current pricing information, go to www.invacare.com/homecare.

Invacare® 3G Arrow® Power Wheelchair w/TrueTrack

BACK SELECTION

ADJASBA BACK TYPE Clear

- SECT Sectional Rear Uprights No Charge
- L8661 Omit Rear Uprights Less \$20.00
- SRUP Straight Rear Uprights No Charge
- L8661 Omit Rear Uprights Less \$20.00
- 9720P No Push Handles (Straight Canes) \$82.00

BACK HEIGHT ADJASBA Clear

- BH12 Back Height 12" No Charge
- BH13 Back Height 13" No Charge
- BH14 Back Height 14" No Charge
- BH15 Back Height 15" No Charge
- BH16 Back Height 16" No Charge
- BH17 Back Height 17" No Charge
- BH18 Back Height 18" No Charge
- BH19 Back Height 19" No Charge
- BH20 Back Height 20" No Charge
- BH21 Back Height 21" No Charge
- BH22 Back Height 22" No Charge
- BH23 Back Height 23" No Charge
- BH24 Back Height 24" No Charge

RECLINER BACK HEIGHT Clear

Used only with manual reclining seat selection.

- BHR18 Back Height 18" No Charge
- BHR20 Back Height 20" No Charge
- BHR22 Back Height 22" No Charge
- BHR24 Back Height 24" No Charge
- BHR26 Back Height 26" No Charge

BACK HEIGHT CONTOURA Clear

- TSBH22 Back Height 22" No Charge
- TSBH25 Back Height 25" No Charge
- TSBH28 Back Height 28" No Charge

BACK ANGLE DEGREE

- 85BA 85° No Charge
- 90BA 90° No Charge
- 95BA 95° No Charge
- 100BA 100° No Charge
- 105BA 105° No Charge

BACK ACCESSORIES Clear

- L8601 Omit Back Upholstery No Charge
- RBP Reinforced Back Upholstery No Charge

ARM SELECTION

ARM TYPE (PAIR)

All arms are flip-up and detachable. Seat widths 12" - 15" not available with fixed height arms must use adjustable height arm style only.

Right Clear

- 18F Full Length Fixed Hgt. Conv. Flip Up No Charge
- 28F Desk Length Fixed Hgt. Conv. Flip Up Standard

Left Clear

- 19F Full Length Adj Hgt. Conv. Flip Up \$230.00 E0973
- 29F Desk Length Adj Hgt. Conv. Flip Up \$230.00 E0973

CANTILEVER ARM TYPE (PAIR) Clear

Cantilever Arms only available with 300 lbs. rated systems. Not available with BPO, Otto Bock, or 12" - 15" seat width. Must use hard back. See chart on last page for arm height adjustment range.

- CLVF Full Length Arm \$350.00
- CLVD Desk Length Arm \$350.00

CANTILEVER MIXED ARM TYPE (EACH)

Right Arm Clear

- CLVFR Full Length Arm Right \$175.00
- CLVDR Desk Length Arm Right \$175.00

Left Arm Clear

- CLVFL Full Length Arm Left \$175.00
- CLVDL Desk Length Arm Left \$175.00

ARMPAD OPTIONS (CANTILEVER ARMS) Clear

- WF Waterfall Armpads Standard
- RCL Reclining Armpads No Charge

NOTE: All specifications and dimensions are approximate.
For the most current pricing information, go to www.invacare.com/homecare.

MK6i™ ELECTRONICS

MK6i CONTROLLER

- MK6TT TrueTrack Controller Standard

EXPANDABLE CONTROLLER SYSTEM

- Clear
- PWH Harness req. for Expandable Electronics \$350.00 E2313
Required to interface joystick or MK6i display with expandable controllers.
- CMPJM6U MPJ+ Multiple Drive Joystick with **Color Display** .. \$750.00 E2377
- ASLRDYM6 MK6 Display Only (No Driver Control) \$750.00 E2377
Display only, No Driver Control. Allows choice from order form.
- MCD Mini Color Display \$750.00 E2377
No port for remote On/Off switch. Stereo port only.

FOR YOUR ASL DRIVE SYSTEM NEEDS
PLEASE REFER TO THE
ASL SINGLE INVOICE PRICE LIST FORM # 02-097.

ALTERNATIVE PROPORTIONAL CONTROL SYSTEMS

- Clear
- 1500M6 Prop RIM Head Control with MK6i Display .. \$2,875.00 E2327
- Clear
- SMHD Small Headrest No Charge
- LGHD Large Headrest No Charge
- 1558MM6 Compact Joystick with Display \$1,130.00 E2373

ALTERNATIVE DIGITAL CONTROL SYSTEMS

- Clear
- Must also use either a Multiple Drive Joystick or MK6i Display.*
- SNPM6 Sip-N-Puff Digital Interface Box \$1,500.00 E2325
- 5020M6 Mini Tash Joystick and Digital Interface Box \$1,750.00 E2321

SIP-N-PUFF KIT

- PKG32666 Therafin® Sip-N-Puff Breath Tube Kit \$380.00 E2326

ADDITIONAL CONTROL CHOICES

- 1558M6 Compact Joystick - Less Std Mount \$770.00
- 5020 Mini Tash Joystick \$405.00
Requires SNPM6 Sip-N-Puff Interface Box.
- PACM6 Prop Attendant Control with quick release mount \$850.00
- 1552M Digital Attendant Control \$850.00
Requires SNPM6 Sip-N-Puff Interface Box.

JOYSTICK MOUNTING HARDWARE

- QLAM6 Quad Link Retract. Joystick Mnt CMPJ+ \$279.00 E1028
Includes ARM250 standard. Not available with TRRO
- QLM6 Quad Link Retract. Joystick Mnt CMPJ+ Van...\$279.00 E1028
Includes ARM250 standard. Not available with TRRO
- ARM250 Stealth Height Adj Joystick Mount \$210.00 E1028
- LEFT Left Handed Mounted Joystick No Charge
- RIMHW RIM Hardware For 1500 \$219.00 E1028
- PKG32669 Complete Bib Assembly (Tash/1558) \$265.00
- GATMPJM6 Gatlin Midline for CMPJ \$600.00 E1028
- GAT1812 Gatlin Midline for 1812 \$600.00 E1028
- RCM Rear Cane Mount for Compact Joystick \$395.00

ELECTRONIC ACCESSORIES

- 1813M6 MK6i Programmer w/Pro Memory Card \$475.00
- 1813SD Pro Memory Card (USB Ready) \$90.00
- 1560 "T" Handle Flexible Joystick Ext. \$77.00 E2323
- 1561 Straight Handle Flexible Joystick Ext. \$77.00 E2323
- 1826 Chin Cup \$48.00 E2324
- A24VPS 24V MK6 Auxiliary Power Source \$200.00
- AUX12M6 Auxiliary Module for outputs 1 & 2 \$795.00
- AUX34M6 Auxiliary Module for outputs 3 & 4 \$795.00
- MPY Monoport "Y" Cable \$15.00
For use with CMPJ+ or MK6i Display to provide additional switch site.
- 2PB Dual Push Button \$147.00
Option to replace Mono Port "Y" cable when requiring two switch functions of the mode switch port on expandable controller systems. Not compatible with SPJ+ Joystick options.
- 2SW Dual Toggle \$172.00
Option to replace Mono Port "Y" cable when requiring two switch functions of the mode switch port on expandable controller systems. Not compatible with SPJ+ Joystick options.
- PBMT Push Button Mounting Kit \$252.00
- RFM Wireless (RF) Mouse Only Module \$750.00
RFM includes Transmitter, USB Radio Frequency Mouse Dongle (Receiver), and Stereo Switch with hook and loop mount for mouse clicks. Not compatible with SPJ+.
- IR Infrared Module with Wireless RF Mouse Control \$995.00
Capable of operating up to six remote Infra-Red controlled devices. Also allows wireless (RF) Mouse operation through the driver control. Compatible with RFM module for control of two wireless mice from one chair. Includes two button stereo switch (with hook/loop mount) for mouse clicks. Not compatible with SPJ+ Joysticks.

POWER LEGS CONTROL MODULE

- Clear
- These options are located in the Configurator under the MK6i Electronics Folder. The Category Header inside the folder is labeled Power Seating Accessories*
- 4WSB Four Way Switch Box No Charge
- Clear
- FWT 4 Way Toggle No Charge
- QPB 4 Push Buttons No Charge
- S4WSB Multiple Actuator Interface Box \$2,600.00 E2311
- Clear
- FWT 4 Way Toggle No Charge
- QPB 4 Push Buttons No Charge

Systems

For systems with two or more actuators, NOT operating through driver control, choose 4WSB with either 4-way toggle or 4-push button.

For systems with 2 or more actuators
OPERATING THROUGH DRIVER CONTROL, choose S4WSB.
Uses 4-way toggle or 4-push button as attendant control.



NOTE: All specifications and dimensions are approximate.
For the most current pricing information, go to www.invacare.com/homecare.

Invacare® 3G Arrow® Power Wheelchair w/TrueTrack

MK6i ELECTRONICS Continued

ASL OPTIONAL ACCESSORIES

- EG5BLK Egg Switch-Black \$156.00
- ASL304A Light Wobble Switch \$180.00
- ASL202J Fiber Optic Switch \$732.00
Must add MK6 Auxiliary Power Source.
- ASL208J Adjustable Proximity Switch \$660.00
Must add MK6 Auxiliary Power Source.
- ASL504M6 Remote Emergency Stop Switch \$1,185.00
- ASL611 Adjustable Switch Mounting Bracket \$252.00

FRAME AND UPHOLSTERY COLOR

CONTEMPORARY

- 24P Wet Black No Charge
- 60P Silver Metallic No Charge
- 105P Midnight Blue No Charge
- 115P Black Prism No Charge

TOO HOT

- 30P Sunny Yellow No Charge
- 61P Electric Red No Charge
- 119P Bubblegum Pink No Charge
- 120P Grape Madness No Charge

WAY COOL

- 62P Electric Blue No Charge
- 127P Grasshopper Green No Charge
- 116P Sunkissed Orange No Charge
- 117P Tangerine Red No Charge
- 29P White Snow No Charge

Discontinued colors 71P Silver Vein, 104P Electric Teal, 121P Lolly Pop Blue, 122P Cosmic Blue, and 125P Emerald Green can be order through specials number C-070109031327

UPHOLSTERY COLORS

- U240 Black Nylon Upholstery Standard
- U67 Black Vinyl No Charge
- GV Grey Vinyl (Van Seat) No Charge

CASTER AND WHEEL SELECTIONS

CASTER OPTIONS

- 1259 6" x 2" Semi-Pneumatic Casters No Charge
Only available with U300 lbs. weight capacity.
- 1228-3 8" x 2 1/4" Semi-Pneu Casters Lt Grey Standard
- 1251-3 8" x 2" Pneumatic Casters No Charge
- 1220-3 9" x 2 3/4" Pneumatic Casters No Charge
Not available with short base. Not available with Transport Tiedown Options.
- 1465-3 9" x 2 3/4" Foam Filled Casters No Charge
Not available with short base. Not available with Transport Tiedown Options.

CASTER SUSPENSION

- Clear
- The below selections are not available with 9" casters.**
- 9213-6 Shock Fork For 6" Casters \$119.00 E1016
- 9213-8 Shock Fork For 8" Casters \$119.00 E1016
Not available with short base. Not available with TRBKTS.

WHEEL OPTIONS

- 1430-3 14" x 3" Pneumatic Tires Standard
- 1431-3 14" x 3" Tires w/Foam Filled Inserts No Charge

REAR WHEEL LOCATION

- RWS Rear Wheel Standard Standard

FENDER OPTION

- L8603 Omit Fender Less \$50.00

POSITIONING STRAPS

SEAT POSITIONING STRAPS

- 1311BK Airline Buckle Seat Positioning Strap \$50.00
Not available with 12" - 15" seat widths.
- 1515 Push Button Style Seat Positioning Strap Standard

CHEST POSITIONING STRAPS

- Clear
- 7311BK Airline Buckle Chest Positioning Strap \$100.00 E0960
Not available with 12" - 15" seat widths.
- 9515 Push Button Style Chest Positioning Strap \$100.00 E0960
- 7321BK Hook and Loop Chest Positioning Strap \$100.00 E0960

ADDUCTORS

- Clear
- Hardware is packaged with equipment and cannot be ordered separately.**
- ADD Adductors Small \$200.00 E0956
- ADDHW Hardware for ADD \$400.00 E1028
- ADDL Adductors Large \$200.00 E0956
- ADDHWL Hardware for ADDL \$400.00 E1028

MISCELLANEOUS

MISCELLANEOUS ACCESSORIES

- VANO2 O2 Holder - Van \$150.00 E2208
- SEATO2 Oxygen Holder (ADJASBA Seat Only) \$150.00 E2208
- 1303 Wheel Locks \$45.00
Not available with pneumatic tires.
- 9205H-3 Heavy Duty Suspension Shock \$400.00

BATTERY SELECTION

BATTERY BOXES

- GP24 Group 24 Battery Boxes Standard

BATTERIES

- 24BATTERY 24 Gel Battery (On Chair) \$410.00 E2363

BATTERY CHARGER OPTION

- Clear
- 110CHARGER 110 Volt Battery Charger No Charge
- LESTR Lester 8 AMP Battery Charger \$320.00 E2367
- L8900 Omit Battery Charger Less \$80.00

**NOTE: All specifications and dimensions are approximate.
For the most current pricing information, go to www.invacare.com/homecare.**

Invacare® 3G Arrow® Power Wheelchair w/TrueTrack

HANGER BRACKET SELECTION

HANGER BRACKET TYPE

Clear

- ASB Center Pivot Style Bracket No Charge
 PSB Pin Style Bracket No Charge
 NHB No Hanger Bracket No Charge

FRONT RIGGING OPTIONS

CENTER PIVOT

Clear

- 60HD 60° Swingaway Footrests 13 1/2" - 17 1/2" .. No Charge
 70NHD 70° Swingaway Footrests 13 1/4" - 17" No Charge
 90SWING 90° Swingaway Footrests No Charge
 70STAPER 70° Tapered Footrests 13" - 17" No Charge
Standard with Impact guards.
 AT5544 Elevating Legrests \$295.00 E0990

PIN STYLE

Clear

- PAS4A Smartleg® Art Legrests 17 1/2" - 21 1/2" \$320.00 K0053
 PW93 Swingaway Footrest 13 1/2" - 19" No Charge
 P904A Swingaway Elevating Legrests 16" - 19 3/4" ... \$255.00 E0990
 PAL4A Smartleg Art Legrests 17 1/4" - 21 1/2" \$320.00 K0053
For seat widths greater than or equal to 16".

POWER LEGS

- ELRPW Power ELR with Articulation 20" (pair) \$1,895.00 E1010
Must pick 4WSB or S4WSB to operate power legs. Independent powered elevating articulating legrests with adjustable angle footplates. Extension range is 13" - 20". Sold as a pair.

FOOTPLATE OPTIONS

Clear

- 1651 Flip-Up Composite Footplates No Charge
 1350 Extra Large Aluminum Footplates No Charge
 1450 Aluminum Footplates-Deluxe No Charge
 AT5543 Adjustable Angle Flip Up Footplates \$165.00 K0040
 AT5043 Flip Up Footplates No Charge

FOOTREST/LEGREST ACCESSORIES

- LPT2 Longer Pivot Slide Tube For 1350 \$83.00
 HLOOP Heel Loops (pair) \$45.00 E0951
 1337 Calf Strap \$27.00 K0038
 IMPACTG Impact Guards \$51.00
 ALPT Longer Pivot and Slide Tube \$84.00
 A9342 "H" Heel Strap \$51.00 K0039

Contoura Seat Sizes							
	D16	D17	D18	D19	D20	D21	D22
W16	X	X	X	X	X	*	*
W18	X	X	X	X	X	*	*
W20	X	X	X	X	X	X	X
W22	*	*	X	X	X	X	X
W24	*		*	X	X	X	X

*Not available for Contoura Seat Cushion.

Back Heights	22", 25" or 28"
--------------	-----------------

Cantilever Arm Height Adjustment			
Back Type	Height	Adjustment	Adjustment with laterals
CURV	20" or greater	11" to 16"	11" to 14"
CURV	19"	11" to 15"	11" to 13"
CURV	18"	11" to 14"	11" to 12"
CURV	17"	11" to 13"	no adjustment
CURV	16"	11" to 12"	not available
Contoura	22", 25" to 28"	12 1/2" to 14 1/2"	no adjustment

Footrest Extension Range			
	AT5043	AT5543	1651
60HD	N/A	12 1/4" - 16 1/2"	13 1/2" - 17 1/2"
70MFX	6 3/4" - 12 1/2"	4 1/2" - 10"	N/A
70NHD	N/A	10 3/4" - 15 3/4"	11 3/4" - 17"
70STAPER	N/A	12 1/2" - 16 1/2"	13 1/4" - 17 1/4"

Measured from Seat Pan to back of Footplate.

Manual Recline Specifications		
Seat Widths	14" - 24"	1" increments
Seat Depths	16" - 24"	1" increments
	(Adjustable 16" - 19", Fixed 20" - 22")	
Back Heights	24" Standard, 18 1/2", 20", 22", 24", 26"	
Headrest	10" Height	
Recline Range	90° to 170° Approximate	
Overall Back Height (without extension)	Standard Base	Heavy-Duty Base
Low Seat	41 1/2"	
Med Seat	43 1/2"	42 3/4"
High Seat	45 1/2"	44.62"

(Measured with 5° Tilt and 10" for Extension)

NOTE: All specifications and dimensions are approximate.
For the most current pricing information, go to www.invacare.com/homecare.

Invacare® 3G Arrow® Power Wheelchair w/TrueTrack

Arrow Specifications		
Maximum Speed	7 MPH	
Seat Width Range (adjustable in 1" increments)	12" - 24"	
Seat Depth Range (adjustable in 1" increments)	12" - 19" Standard Base	20" - 22" Long Base
Seat-to-Floor	Standard Base	
	17 1/2" Standard	
	19 3/4" Optional	
	21" Optional	
Overall Length	29 1/2"	Short Base
	32 1/2"	Long Base
Overall Width	25"	
Weight Limits	300 lbs.	
	400 lbs.	
Measurement based on wheel base. Overall width is width picked plus 7 1/2" for armrest and joystick. Measured from headtube to anti-tipper. Based on 5° Tilt and 16" Seat Depth. Add 1/4" for additional 2" of Seat Depth. Field adjustable for 10° of total pre-tilt. Note: All specifications and dominions are approximate.		

Chair Features and Options

- Independent Suspension for RWD for Clinical and Functional Comfort
- Full Electronic Communication & Switch Option Capability
- Programmable Performance to Meet Individuals Needs
- Four Drives 1, 2, 3 & 4 - Four Chairs in One
- 5" Ground Clearance (RWD)
- Solid Seat Pan
- Convenient Rear Battery Access
- Automatic Self Connecting Battery Hookup
- Adjustable Depth Rigging Attachment 2" {1" increments}
- Automatic Dual Mode Battery Charger (Lester Brand)
- Multiple Colors to Choose From
- Option of Adjustable Tilt and Back Angle Seat

Weight Limit 400 Pounds

Invacare recommends all wheelchairs be fitted and programmed to the consumer under the guidance of a healthcare professional, knowledgeable of the particular limitations of the product.

Note: All weight limits are inclusive of add-on equipment and seat systems.

Memory Card - BASIC

- Standard with ALL Expandable Control Systems (CMPJ+, PSR+, PSF+, MK6i Display).
- Used to back up / restore programmed settings for only one chair.
- May be given to user or kept with provider's files for safe keeping.
- Does not contain advanced diagnostics / help screens / file structure.
- Not compatible with SPJ+ Joysticks.

Memory Card - PROFESSIONAL

- Standard with ALL MK6i Programmers.
- Available also with USB Card Reader on Order Form.
- Contains advanced diagnostics, help screens, software upgrades, file storage / retrieval
- Used to back up multiple chairs and programming settings.
- Intended for use by qualified providers only.
- Not compatible with SPJ+ Joysticks.

Seating and Positioning

For your seating needs, please refer to the following rehab seating and positioning products: Price List Form No. 07-027

Comments

NOTE: All specifications and dimensions are approximate.
For the most current pricing information, go to www.invacare.com/homecare.



Invacare Corporation

USA

One Invacare Way

Elyria, OH USA

44035-4190

800-333-6900

www.invacare.com

©2013 Invacare Corporation

All rights reserved. Trademarks are identified by the symbols™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications and prices are subject to change without notification.