Invacare *e*Forms

For additional information, visit www.invacare.com/eForms

Customer Service: 800-333-6900

Fax: 800-678-4682

www.invacare.com To order Literature: 800-828-6282



Invacare *e*Forms offer new functionality!

Invacare is pleased to offer you an upgraded solution to your ordering process. Our enhanced order forms allow you to fill out a form electronically, print and fax the form, save and email* it to Customer Service, or maintain the business practices that work for you today. The format has been revised to reveal a cleaner look with electronic selection and input functions. *Saving and emailing the eForms is available with Acrobat Professional software

If you are using Adobe Reader

- View and download all forms
- Save blank versions of the forms to your desktop
- Complete forms electronically
- Print blank or completed forms
- Fax completed forms to Customer Service

Please note: The 'Email' button at the end of the form will not work with Adobe Reader. The 'Email' button will only work while using Adobe Acrobat Professional.

If you are using Adobe Acrobat Professional

- View and download all forms
- **Save** blank AND completed versions of the forms to maintain electronic records
- Complete forms electronically
- Print blank or completed forms
- Fax completed forms to Customer Service
- Email completed forms to Customer Service by using the 'Email' button at the end of the form.

For information on downloading Adobe Reader for free, or puchasing Adobe Acrobat Professional, please visit <u>www.adobe.com</u>.

Invacare[®] Formula[™] CG Tilt Only for Storm Series[®] Power Wheelchair

Price List and Order Form Customer Service: 800-333-6900	Fax: 800-678-4682	To order	www.invacare.c literature: 800-828-628 YPE: Quote	32	Yes, you can:
Date Account #			P.O. #		
Company Name					
Address					
Contact Name		Comments			
		Special No	te		
HCPCS codes are not intended to be, nor submitting claims to the Medicare Program		ng or legal advice. Pr	oviders are responsible		
POWERED SEATING BASE [™] SP for Formula CG with MK6i [™] \$6,850.00 K0877 The Torque SP Base comes Standard with High Performance 4-Pole Motors, and up to a top speed of 6 1/2 MPH. 3GTQ3-CG Torque 3 for Formula CG with MK6i \$5,970.00 The Torque 3 Base comes Standard with 4-Pole Motors and up to a top speed of 5.8 MPH.		0.00 ^{K0856}	TRANSPORT T	Seat Frame for Base IE DOWN Ichair Transport Bra	e Only Less \$600.00 ickets \$250.00 ALL seat-to-floor height. Not available
Group 4 - K0877 B	asic Package			TRANSPORT T	IE DOWN*
3GTQ-CG Torque SP for Formula CG with MK6i 300 lbs. Weight Capacity 1431-3 14" x 3" Tires with Foam Filled Inserts No Charge Performance 4 Pole Motors Standard 1228-3 8" x 2 1/4" Semi Pneumatic Casters Standard Flip Back Adjustable Height Arm Standard MK690ACC Controller with Actuator Control Standard 70NHD 70° Swingaway Footrests Standard 8 Amp Off Board Battery Charger Standard 1515 Push Button Style Seat Positioning Strap Standard Seat Option CGT Tilt (E1002) \$4,590.00		lard lard lard lard lard lard lard	TRBKTS has not b transport brackets on As of this date, the tie-down systems f moving vehicle of a wheelchairs should for transportation a the auto industry. I wheelchair transpo	four factory installe een crash-tested in a ly to secure an unoccu Department of Tran or transportation of a any type. It is Invaca be transferred into a nd use be made of the nvacare cannot and of rtation systems.	d wheelchair transport brackets. accordance with WC19. Use these pied wheelchair during transport. Isportation has not approved any a user while in a wheelchair, in a are's position that users of appropriate seating in vehicles he restraints made available by does not recommend any
Group 3 - K0856 B	asic Package		FRAME TYPE		

Group 3 - K0856 Basic Packag	e
3GTQ3-CG Torque 3 for Formula CG with MK6i	
300 lbs. Weight Capacity	
1431-3 14" x 3" Tires with Foam Filled Inserts	No Charge
4 Pole Motors	Standard
1228-3 8" x 2 1/4" Semi Pneumatic Casters	Standard
Flip Back Adjustable Height Arm	Standard
MK690ACC Controller with Actuator Control	Standard
70NHD 70° Swingaway Footrests	Standard
8 Amp Off Board Battery Charger	Standard
1515 Push Button Style Seat Positioning Strap	Standard
Seat Option	
CGT Tilt (E1002)	\$4,590.00

SFS Short Base No Charge Not Available on seat depths of 1922D. Must pick Long Base. Long Base No Charge

CASTER OPTIONS

1228-3	8" x 2 1/4" Semi-Pneu Casters Lt Grey	Standard
○ 1251-3	8" x 2" Pneumatic Casters	No Charge

CASTER SUSPENSION

09213-8	Shock Fork For 8" Casters		E1016
	Not available with short base.	Not available with TRBKTS.	

DRIVE WHEEL OPTIONS

1430-3	14" x 3" Pneumatic Tires	Standard
1431-3	14" x 3" Tires with Foam Filled Inserts	No Charge

FRAME FINISH

CONTEMPORARY

○ 24P	Wet Black	No Charge
○ 60P	Silver Metallic	No Charge
	Midnight Blue	
○ 115P	Black Prism	No Charge

TOO HOT

○ 30P	Sunny Yellow	No Charge
○ 61P	Electric Red	No Charge
○ 119P	Bubblegum Pink	No Charge
C 120P	Grape Madness	No Charge

WAY COOL

○ 62P	Electric Blue	No Charge
○ 127P	Grasshopper Green	No Charge
○116P	Sunkissed Orange	No Charge
○ 117P	Tangerine Red	No Charge
C 29P	White Snow	No Charge

SEAT

SYSTEM TYPE

○ CGT	CG Tilt Only System	\$4,590.00 E1002
\bigcirc	NVACT Conventional Actuator	Standard

SEAT SELECTION

○ CTB	Contoura [™] Seating \$605.00	E2620 E2621
O BPO	Does not include seat cushion. Back Post Only, Round Canes Only No Charge	
	Push handles included.	
○ CURV	Curved No Charge	
	Push handles included.	

ADJUSTABLE SEAT WIDTH RANGE

○ 1620W	16" to 20" Seat Width (Adjustable) No Charge	
○ 2024W	20" to 22" Seat Width (Adjustable) No Charge	

SEAT WIDTH SETTING

○W16	16" wide	No Charge
○W17	17" wide	No Charge
W18	18" wide	No Charge
∕W19	19" wide	No Charge
○ W20	20" wide	No Charge
○W21	21" wide	No Charge
○W22	22" wide	No Charge

SEAT-TO-FLOOR HEIGHT

Due to variation	in options the listed seat-to-floor heights have	e a tolerance of +/5"
CLOW	$18.50"$ at 0° Not available with Transport Tie Down opt	

○ MED	19.50" at 0° No C	Charge
○ TALL	20.50" at 0° No 0	Charge

ADJUSTABLE SEAT DEPTH RANGE

🔿 1619D	16" to 19" Seat Depth Range	. No Charge
() 1922D	19" to 22" Seat Depth Range	No Charge

ADUSTABLE SEAT DEPTH SETTING

OD16	Seat Depth 16"	No Charge
OD17	Seat Depth 17"	No Charge
O D18	Seat Depth 18"	No Charge
O D19	Seat Depth 19"	No Charge
O D20	Seat Depth 20"	No Charge
OD21	Seat Depth 21"	No Charge
OD22	Seat Depth 22"	No Charge

SEAT TILT

Indicates static tilt of seat. Impacts seat-to-floor height. Every 5° raises the seat-tofloor boight

noor neight.		
○ 0SA	0° Sta	andard

BACK HEIGHT

◯ TSBH16	Back Height 16"	No Charge
◯ TSBH17	Back Height 17"	No Charge
◯ TSBH18	Back Height 18"	No Charge
◯ TSBH19	Back Height 19"	No Charge
◯ TSBH20	Back Height 20"	No Charge
◯ TSBH21	Back Height 21"	No Charge
◯ TSBH22	Back Height 22"	No Charge
◯ TSBH23	Back Height 23"	No Charge
◯ TSBH24	Back Height 24"	No Charge
◯ TSBH25	Back Height 25"	No Charge
◯ TSBH26	Back Height 26"	No Charge

CONTOURA BACK HEIGHT

◯ TSBH22	Back Height 22"	No Charge
◯ TSBH25	Back Height 25"	No Charge
TSBH28	Back Height 28"	No Charge

			Short Frame (16-19" Seat Depth)			Long Frame (19-22" Seat Depth)		
			Low	Medium	Tall	Low	Medium	Tall
Powered	Tilt, Tilt/	0 Degree						
Seating	Recline	Seat Angle	18.5	19.5	20.5	18.75	19.75	20.75

NOTE: All specifications and dimensions are approximate.

For the most current pricing information, go to www.invacare.com/homecare.

SEAT Continued

ADJUSTABLE ARM TYPE (PAIR)

🔿 19F	Full Length Adj. Hgt	No Charge
○ 29F	Desk Length Adj. Hgt	No Charge

MIXED ADJUSTABLE ARM TYPE (EACH)

Right Arm		
19FR	Full Length Adj. Hgt. Right	. No Charge
29FR	Desk Length Adj. Hgt. Right	. No Charge
Left Arm		
19FL	Full Length Adj. Hgt. Left	. No Charge
29FL	Desk Length Adj. Hgt. Left	. No Charge

SEAT POSITIONING STRAPS

○ 1515	Push Button Style Seat Positioning Strap Standard	
() 1515	Push Button Style Seat Positioning Strap Standard	

CHEST POSITIONING STRAPS

POSITIONING ACCESSORIES

HEADREST TYPE (OPTIONAL)

Hardware is packaged with equipment and cannot be ordered separately.

C TAR127700 Curved Headrest \$210.00	E0955
TAR127700HW Removable & Adjustable Hardware \$210.00	E1028

ADDUCTORS

Hardware is packaged with equipment and cannot be ordered separately.

🔿 ADD	Addu	ctors Small	 \$200.00	E0956
\bigcirc	ADDHW	Hardware for ADD	 \$400.00	E1028

LATERALS (Pair)

Curved Back / BPO Hardware is packaged with equipment and cannot be ordered separately.

TCSLM Clamp-On Swingaway Laterals Medium \$210.00 EC	0956
○ TCSLMHW Removable and Adjustable Hardware \$210.00 E	1028
TCSLL Clamp-On Swingaway Laterals Large \$210.00 EC	0956
○ TCSLLHW Removable and Adjustable Hardware \$210.00 FI	1028

Contoura Back

For use with Contoura seating only. Hardware is packaged with equipment and cannot be ordered senar:

and bannet be endered coparatory.	
CSLM Swingaway Contoura Lat. Med \$210.00 E	E0956
CSLMHW Removable and Adjustable Hardware \$210.00 E	E1028
CSLL Swingaway Contoura Lat. Lg \$210.00 E	E0956
CSLLHW Removable and Adjustable Hardware \$210.00 E	E1028

ELECTRONICS

MK6i CONTROLLERS Torque SP / Torque 3

MK690ACC 90 AMP Controller with Actuator Contol Standard

G-Trac Module \$995.00 K0108 GTRAC G-Track provides Improved Tracking Capabilities and Enhanced Traction. Does not require Special Motors. Not available with SPJ+ joysticks or MK6TT.

STANDARD DRIVER CONTROL

Torque SP / Torque 3

\odot	SPJAP	SPJ	+ with Power	Seating	Switch	Standard
---------	-------	-----	--------------	---------	--------	----------

ACTUATOR CONTROL THROUGH SPJ+

○ SPJAPM SPJ+ with Joystick Actuator Control \$1,300.00 E2310 Optional for Torque SP and Torque 3. Allows actuator operation through Driver Control. Includes port for attendant . switch.

EXPANDABLE CONTROLLER SYSTEM

○ PWH	Harness req. for Expandable Electronics \$350.00 E2313
	Required to interface joystick or MK6i display with expandable
	controllers.

C ASLRDYM6 MK6i Display Only (No Driver Control) \$750.00 E2377 Display only, No Driver Control. Allows choice from order form.

ALTERNATIVE PROPORTIONAL CONTROL SYSTEMS

○ 1558MM6 Compact Joystick with Display \$1,130.00 E2373

ALTERNATIVE DIGITAL CONTROL SYSTEMS Must also use either a Multiple Drive Joystick or MK6i Display

○ SNPM6	Sip-N-Puff Digital Interface Box \$1,500.00 E2325	5

SIP-N-PUFF KIT

ADDITIONAL CONTROL CHOICES

\bigcirc	1558M6	Compact Joystick - Less Std Mount	
------------	--------	-----------------------------------	--

PACM6 Prop Attendant Control with quick release mount \$850.00

JOYSTICK MOUNTING HARDWARE

O QLAM6	Quad Link Retract. Joystick Mnt CMPJ+ \$279.00 Includes ARM20 Standard.	E1028
O QLA	Quad Link Retract. Joystick Mount SPJ+ \$279.00 Includes ARM20 Standard.	E1028
~	Stealth Height Adjustable Joystick Mount \$210.00	E1028
\sim	Left Handed Mounted Joystick No Charge	
ELECTRO	NIC ACCESSORIES	
O 1813M6	MK6i Programmer w/Pro Memory Card \$475.00	
1912CD	Pro Mamory Card (USPPaady) \$00.00	

I

🔿 1813M6	MK6i Programmer w/Pro Memory Card	\$475.00
🔿 1813SD	Pro Memory Card (USBReady)	\$90.00
O AUX12M6	Auxiliary Module for outputs 1 & 2	\$795.00
O AUX34M6	Auxiliary Module for outputs 3 & 4	\$795.00
○ MPY	Monoport "Y" Cable for CMPJ+ or Display.	\$15.00
○ 2PB	Dual Push Button	\$147.00
	Option to replace Mono Port "Y" cable when requir switch functions of the mode switch port on expan controller systems. Not compatible with SPJ+ Joys	dable
O PBMT	Push Button Mounting Kit	\$252.00
C A24VPS	MK6i Auxiliary Power Source	\$200.00

POWERED SEATING ACCESSORIES

SANODE One Actuator Operation thru Driver Cntrl \$1,300.00 E2310 Required for Tilt Operation through Driver Control - Single Actuator Systems Only. Not available with SPJ+ Joystick. Does not come with Egg switch.			
○ 4WSB Four	Way Switch Box	No Charge	
○ FWT	4 Way Toggle	No Charge	
⊂ QPB	4 Push Buttons	No Charge	
S4WSB Mult	iple Actuator Interface Box		
○ FWT	4 Way Toggle	No Charge	
⊂ QPB	4 Push Buttons	No Charge	

Systems

For systems with two or more actuators, NOT operating through driver control, choose 4WSB with either 4-way toggle or 4-push button.

For systems with 2 or more actuators OPERATING THROUGH DRIVER CONTROL, choose S4WSB. Uses 4-way toggle or 4-push button as attendant control.

NOTE: All specifications and dimensions are approximate. For the most current pricing information, go to www.invacare.com/homecare.

ELECTRONICS Continued

ASL OPTIONAL ACCESSORIES

C EGSBLK	Egg Reset-Black	\$156.00
ASL611	Adjustable Switch Mounting Bracket	\$252.00

FRONT RIGGINGS

HANGER BRACKET TYPE

CASB Center Pivot Style Bracket No Charge

CENTER PIVOT STYLE (PAIR)

○ 70NHD	70° Swingaway Footrests 13 1/4"-17" No Charge	
○ 70STAPER	70° Tapered Footrests 13"-17" No Charge	
	Not available with 19" or 20" wide and 9"casters. Standard with	
_	impact guards on the footrests.	
○ AT5044	Elevating Legrests \$255.00	E0990
\mathbf{O}	Elevating Legrests	

FOOTPLATE OPTIONS (PAIR)

1651	Flip-Up Composite Footplates No Charge	
1350	Extra Large Aluminum Footplates No Charge	
O AT5543	Adjustable Angle Flip-Up Footplates \$165.00	K0040

FOOTREST/LEGREST ACCESSORIES

1337	Calf Strap	\$27.00	E0038
○ IMPACTG	Impact Guards	\$51.00	
○ HLOOP	Heel Loops (pair)	\$45.00	E0951
C LPT2	Longer Pivot Slide Tube For 1350	\$83.00	
○ ALPT	Longer Pivot And Slide Tube	\$84.00	
C ALPT4	Longer Pivot And Slide Tube 4"	\$84.00	

MISCELLANEOUS

BATTERY BOXES

 GP24
 Group 24 Battery Boxes
 No Charge

 Must select TRBKTS.
 Must select TRBKTS.
 Must select TRBKTS.

BATTERIES

Batteries will be installed by Invacare Corporation at no additional charge. Chair	
will be shipped complete with batteries.	

○ 24BATTERY	24 Gel Battery (On Chair)	\$410.00 E2363
-------------	---------------------------	----------------

BATTERY CHARGER OPTION

110CHARGE	R 110 Volt Battery Charger	No Charge
C L8900	Omit Battery Charger	Less \$80.00

MISCELLANEOUS ACCESSORIES

1303	Wheel Locks	\$45.00
09205-3	Suspension Shocks	\$400.00 E1016
◯ FDR	Fender Not available with chair with TRBKTS.	\$320.00

	Contoura Seat Sizes						
	D16	D17	D18	D19	D20	D21	D22
W16	*	*	*	*	*	*	*
W18	*	*	*	*	*	*	*
W20	*	*	*	*	*	*	*
W22	*	*	*	*	*	*	*
*Not ava	*Not available with Contoura Seat Cushion						
Back Heights 22", 25" or 28"							
Curved							
Seat Widths			16"- 24"				
Seat Depths			16"- 22"				
Back Heights			16"- 26"				

Back Height Chart			
Back Type Back Height			
CTB - Contoura	22", 25" and 28"		
BPO - Back Posts Only	20"-24"		
CURV - Curved	16"-26"		

Battery Boxes			
Standard Battery Tray*			
Torque SP/Torque 3 Group 24			
*Sealed Gel Cell Batteries Recommended			

Memory Card - BASIC

• Standard with ALL Expandable Control Systems (CMPJ+, MK6i Display).

- Used to Back Up/Restore Programmed Settings for Only One Chair.
- May be given to user or kept with provider's files for safe keeping.
- Does not Contain Advanced Diagnostics/Help Screens/File Structure.
- Not compatible with SPJ+ Joysticks.

Memory Card - PROFESSIONAL

- Standard with ALL MK6i Programmers.
- Available also with USB Card Reader on Order Form.
- Contains Advanced Diagnostics, Help Screens, Software Upgrades,
- File Storage/Retrieval

Yes, you can. Revised December 1, 2016

- Used to back up multiple chairs and programming settings.
- Intended for use by qualified providers only.
- Not compatible with SPJ+ Joysticks.

Batteries

GP24 batteries with terminal configuration (negative on the left and positive on the right) MUST be used. GP24 batteries that have the reverse terminal configuration MUST not be used. Terminals MUST have a cross hole for proper battery connection. See Owner's Manual, part number 1114809. These recommendations MUST be followed otherwise injury and damage may occur.

Tech Tips

Operating Powered Seating Through the Driver's Control. All chairs require adding or replacing one component to allow actuator operation through the driver control.

SPJ+ with PSS: (SPJAP)

Replace joystick with SPJ+ with ACC (Actuator Control) (SPJAPM) - E2310

Single Actuator Systems with CMPJ+, MK6i Display: Add a SANODE (Single Actuator Node) to the system - E2310

Multiple Actuator Systems:

Replace the standard 4 Way Switch box with a Multiple Actuator Interface Box - E2311

Formula CG Tilt Only Features and Specifications

- Weight capacity of 300 lbs. or 400 lbs. (with heavy-duty motors).
- Curved, Conventional or Contoura style seating.
- Adjustable seat frame with Conventional style.
- Seat-to-Floor height of 17 3/4" at 0° tilt.

Seating and Positioning

For your seating needs, please refer to the following rehab seating and positioning products: Price List Form No. 07-027

Comments

NOTE: All specifications and dimensions are approximate. For the most current pricing information, go to www.invacare.com/homecare.

Invacare Corporation USA

One Invacare Way Elyria, OH USA 44035-4190 800-333-6900

www.invacare.com

©2016 Invacare Corporation All rights reserved. Trademarks are identified by the symbols[™] and [®]. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications and prices are subject to change without notification.