



Yes, you can.®

Invacare eForms offer new functionality!

Invacare is pleased to offer you an upgraded solution to your ordering process. Our enhanced order forms allow you to fill out a form electronically, print and fax the form, save and email* it to Customer Service, or maintain the business practices that work for you today. The format has been revised to reveal a cleaner look with electronic selection and input functions.

*Saving and emailing the eForms is available with Acrobat Professional software

If you are using Adobe Reader

- **View** and download all forms
- **Save** blank versions of the forms to your desktop
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service

Please note: The 'Email' button at the end of the form will not work with Adobe Reader. The 'Email' button will only work while using Adobe Acrobat Professional.

If you are using Adobe Acrobat Professional

- **View** and download all forms
- **Save** blank AND completed versions of the forms to maintain electronic records
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service
- **Email** completed forms to Customer Service by using the 'Email' button at the end of the form.

For information on downloading Adobe Reader for free, or purchasing Adobe Acrobat Professional, please visit www.adobe.com.

Please consider the environment before printing this page with your order form.
For more information on features and benefits of electronic order forms, please visit www.invacare.com/eForms

Invacare® Formula™ CG Tilt/Recline for Storm Series® Power Wheelchair

Price List and Order Form

www.invacare.com

Customer Service: 800-333-6900

Fax: 800-678-4682

To order literature: 800-828-6282



Yes, you can.®

TRANSACTION TYPE:

- Order Quote

Date _____ Account # _____ P.O. # _____
 Company Name _____ Phone # _____
 Address _____ City _____ State _____ Zip Code _____ Country _____
 Contact Name _____ Comments _____

Special Note

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.

POWERED SEATING BASE

- 3GTQ-MCG Torque™ SP for Formula CG with MK6i™ .. \$7,700.00 ^{K0884}
The Torque SP base come standard with performance 4-pole motors and a top speed of 6 1/2 mph.
- 3GTQ3-MCG Torque 3 for Formula CG with MK6i™ \$6,275.00 ^{K0861}
The Torque 3 comes standard with 4-pole motors and a top speed of 5.8 mph.

OMIT SEAT FRAME

- OSF Omit Seat Frame for Base Only Less \$600.00

TRANSPORT TIE DOWN

- TRBKTS Wheelchair Transport Brackets \$250.00
Only available with MED or TALL seat-to-floor height. Not available with vent tray.

TRANSPORT TIE-DOWN*

TRBKTS includes four factory installed wheelchair transport brackets. TRBKTS has not been crash-tested in accordance with WC19. Use these transport brackets only to secure an unoccupied wheelchair during transport. As of this date, the Department of Transportation has not approved any tie-down systems for transportation of a user while in a wheelchair, in a moving vehicle of any type. It is Invacare's position that users of wheelchairs should be transferred into appropriate seating in vehicles for transportation and use be made of the restraints made available by the auto industry. Invacare cannot and does not recommend any wheelchair transportation systems.

* Transport Tie Down option is not retrofittable to existing models and is not field serviceable.

FRAME TYPE

- SFS Short Base No Charge
Not available with seat depths 19"- 22". Must pick long base.
- LFS Long Base No Charge

Group 4 - K0884 Basic Package

3GTQ-MCG Torque SP for Formula CG with MK6i	
300 lbs. Weight Capacity	
1431-3 14" X 3" Tires with Foam Filled Inserts	No Charge
Performance 4-Pole Motors	Standard
1228-3 8" X 2 1/4" Semi Pneumatic Casters	Standard
Flip Back Adjustable Height Arm	Standard
MK690ACC Controller with Actuator Control	Standard
70NHD 70° Swingaway Footrests	Standard
8 Amp Off Board Battery Charger	Standard
1515 Push Button Style Seat Positioning Strap	Standard
Seat Option	
CGTR Tilt/Recline (E1007)	\$9,595.00

Group 3 - K0861 Basic Package

3GTQ3-MCG Torque 3 for Formula CG with MK6i	
300 lbs. Weight Capacity	
1431-3 14" x 3" Tires with Foam Filled Inserts	No Charge
4 Pole Motors	Standard
1228-3 8" x 2 1/4" Semi Pneumatic Casters	Standard
Flip Back Adjustable Height Arm	Standard
MK690ACC Controller with Actuator Control	Standard
70NHD 70° Swingaway Footrests	Standard
8 Amp Off Board Battery Charger	Standard
1515 Push Button Style Seat Positioning Strap	Standard
Seat Option	
CGTR Tilt/ Recline (E1007)	\$9,595.00

CASTER OPTIONS

- 1228-3 8" x 2 1/4" Semi-Pneu Casters Lt Grey Standard
- 1251-3 8" x 2" Pneumatic Casters No Charge

CASTER SUSPENSION

- 9213-8 Shock Fork For 8" Casters \$119.00 ^{E1016}
Not available with short base. Not available with TRBKTS.

WHEEL OPTION

- 1430-3 14" x 3" Pneumatic Tires Standard
- 1431-3 14" x 3" Tires with Foam Filled Inserts No Charge

**NOTE: All specifications and dimensions are approximate.
 For the most current pricing information, go to www.invacare.com/homecare.**

Invacare® Formula™ CG Tilt/Recline for Storm Series® Power Wheelchair

FRAME FINISH

CONTEMPORARY

- 24P Wet Black No Charge
- 60P Silver Metallic No Charge
- 105P Midnight Blue No Charge
- 115P Black Prism No Charge

TOO HOT

- 30P Sunny Yellow No Charge
- 61P Electric Red No Charge
- 119P Bubblegum Pink No Charge
- 120P Grape Madness No Charge

WAY COOL

- 62P Electric Blue No Charge
- 127P Grasshopper Green No Charge
- 116P Sunkissed Orange No Charge
- 117P Tangerine Red No Charge
- 29P White Snow No Charge

SEAT

SYSTEM TYPE

- CGTR Tilt/Recline w/Mechanical Shear Reduction . \$9,595.00 E1007

ACTUATOR

- CNVACTS Conventional Actuators Standard

SEAT SELECTION

- CNB Conventional Back No Charge
- CTB Contoura™ Seating \$605.00 E2620 E2621
Does not include seat cushion

ADJUSTABLE SEAT WIDTH RANGE

- 1620W 16" to 20" Seat Width Range (Adjustable) No Charge
- 2024W 20" to 22" Seat Width Range (Adjustable) No Charge

SEAT WIDTH SETTING

- W16 16" Seat Width No Charge
- W17 17" Seat Width (Not available with Contoura) No Charge
- W18 18" Seat Width No Charge
- W19 19" Seat Width (Not available with Contoura) No Charge
- W20 20" Seat Width No Charge
- W21 21" Seat Width (Not available with Contoura) No Charge
- W22 22" Seat Width No Charge

ADJUSTABLE SEAT DEPTH RANGE

- 1619D 16" to 19" Seat Depth Range (Adjustable) No Charge
- 1922D 19" to 22" Seat Depth Range (Adjustable) No Charge

SEAT DEPTH SETTING

- D16 16" Seat Depth No Charge
- D17 17" Seat Depth No Charge
- D18 18" Seat Depth No Charge
- D19 19" Seat Depth No Charge
- D20 20" Seat Depth No Charge
- D21 21" Seat Depth No Charge
- D22 22" Seat Depth No Charge

SEAT-TO-FLOOR HEIGHT

Due to variation in options the listed seat-to-floor heights have a tolerance of +/- .5"

- LOW 18.50" at 0° No Charge
Not available with Transport Tie Down option.
- MED 19.50" at 0° No Charge
- TALL 20.50" at 0° No Charge

SEAT TILT

Indicates static tilt of seat, impacts seat-to-floor height. Every 5 degrees raises the seat-to-floor height.

- 0SA 0° Standard

CONVENTIONAL BACK HEIGHT

- TSBH20 Back Height 20" No Charge
- TSBH21 Back Height 21" No Charge
- TSBH22 Back Height 22" No Charge
- TSBH23 Back Height 23" No Charge
- TSBH24 Back Height 24" No Charge
- TSBH25 Back Height 25" No Charge
- TSBH26 Back Height 26" No Charge

CONTOURA BACK HEIGHT

- TSBH25 Back Height 25" No Charge
- TSBH28 Back Height 28" No Charge

ADJUSTABLE ARM TYPE (EACH)

Right Arm

- RA19FR Full Length Adj. Hgt. Right No Charge
- RA29FR Desk Length Adj. Hgt. Right No Charge

Left Arm

- RA19FL Full Length Adj. Hgt. Left No Charge
- RA29FL Desk Length Adj. Hgt. Left No Charge

SEAT POSITIONING STRAPS

- 1515 Push Button Style Seat Positioning Strap Standard

CHEST POSITIONING STRAPS

- 7321BK Hook & Loop Chest Positioning Strap \$100.00 E0960

			Short Frame (12-16" or 16-19" Seat Depth)			Long Frame (19-22" Seat Depth)		
			Low	Medium	Tall	Low	Medium	Tall
Powered Seating	Tilt, Tilt/Recline	0 Degree Seat Angle	18.5	19.5	20.5	18.75	19.75	20.75
		5 Degree Seat Angle	20.75	21.75	22.75	21.25	22.25	23.25

NOTE: All specifications and dimensions are approximate.
For the most current pricing information, go to www.invacare.com/homecare.

Invacare® Formula™ CG Tilt/Recline for Storm Series® Power Wheelchair

POSITIONING ACCESSORIES

HEADREST TYPE (OPTIONAL)

Hardware is packaged with equipment and cannot be ordered separately.

- TARI27700 Curved Headrest \$210.00 E0955
- TARI27700HW Removable & Adjustable Hardware \$210.00 E1028

LATERALS (Pair)

CONVENTIONAL BACK

Hardware is packaged with equipment and cannot be ordered separately.

- TMSLM Swingaway Laterals Medium \$210.00 E0956
 - TMSLMHW Removable and Adjustable Hardware \$210.00 E1028
- TMSLL Swingaway Laterals Large \$210.00 E0956
 - TMSLLHW Removable and Adjustable Hardware \$210.00 E1028

CONTOURA SEATING

For use with Contoura seating only. Hardware is packaged with equipment and cannot be ordered separately.

- CSLM Swingaway Contoura Lat. Med \$210.00 E0956
 - CSLMHW Removable and Adjustable Hardware \$210.00 E1028
- CSLL Swingaway Contoura Lat. Lg \$210.00 E0956
 - CSLLHW Removable and Adjustable Hardware \$210.00 E1028

ADDUCTORS

Hardware is packaged with equipment and cannot be ordered separately.

- ADD Adductors Small \$200.00 E0956
- ADDHW Hardware for ADD \$400.00 E1028

ELECTRONICS

MK6i CONTROLLERS

Torque SP / Torque 3

- MK690ACC 90 AMP Controller with Actuator Control Standard
- GTRAC G-Trac Module \$995.00 K0108
G-Trac provides Improved Tracking Capabilities and Enhanced Traction. Does not require Special Motors. Not available with SPJ+ joysticks.

EXPANDABLE CONTROLLER SYSTEM

- PWH Harness req. for Expandable Electronics \$350.00 E2313
Required to interface joystick or MK6i display with expandable controllers.
- CMPJM6 MPJ™+ Multiple Drive Joystick with Color Display \$750.00 E2377
- ASLRDYM6 MK6i Display Only (No Driver Control) \$750.00 E2377
Display only, No Driver Control. Allows choice from order form.

FOR INFORMATION ON ASL ELECTRONICS CALL
1-800-626-8698 OR VISIT WWW.ASL-INC.COM

ELECTRONICS Continued

ALTERNATIVE PROPORTIONAL CONTROL SYSTEMS

- 1558MM6 Compact Joystick with Display \$1,130.00 E2373

ALTERNATIVE DIGITAL CONTROL SYSTEMS

Must also use either a Multiple Drive Joystick or MK6 Display.

- SNPM6 Sip-N-Puff Digital Interface Box \$1,500.00 E2325

SIP-N-PUFF KIT

- PKG32666 Therafin Sip-N-Puff Breath Tube Kit \$380.00 E2326

ADDITIONAL CONTROL CHOICES

- 1558M6 Compact Joystick - Less Std Mount \$770.00
- PACM6 Prop Attendant Control with quick release mount \$850.00

JOYSTICK MOUNTING HARDWARE

- QLAM6 Quad Link Retract. Joystick Mnt CMPJ+ \$279.00 E1028
Includes ARM250 standard.
- ARM250 Stealth Height Adjustable Joystick Mount \$210.00 E1028
- LEFT Left Handed Mounted Joystick No Charge

ELECTRONIC ACCESSORIES

- 1813M6 MK6i Programmer with Pro Memory Card \$475.00
- 1813SD Pro Memory Card (USB ready) \$90.00
- A24VPS 24V MK6i/Auxiliary Power Source \$200.00
- AUX12M6 Auxiliary Module for outputs 1 & 2 \$795.00
- AUX34M6 Auxiliary Module for outputs 3 & 4 \$795.00
- 2PB Dual Push Button \$147.00
Option to replace Mono Port "Y" Cable when requiring two switch functions of the mode switch port on expandable controller systems. Not compatible with SPJ+ Joystick options.
- PBMT Push Button Mounting Kit \$252.00
- MPY Monoport "Y" Cable for CMPJ+ or Display \$15.00

POWERED SEATING ACCESSORIES

- 4WSB Four Way Switch Box No Charge
 - FWT 4 Way Toggle No Charge
 - QPB 4 Push Buttons No Charge
- S4WSB Multiple Actuator Interface Box \$2,600.00 E2311
 - FWT 4 Way Toggle No Charge
 - QPB 4 Push Buttons No Charge

Systems

For systems with two or more actuators, NOT operating through driver control, choose 4WSB with either 4-way toggle or 4-push button.

For systems with 2 or more actuators OPERATING THROUGH DRIVER CONTROL, choose S4WSB.
Uses 4-way toggle or 4-push button as attendant control.

NOTE: All specifications and dimensions are approximate.
For the most current pricing information, go to www.invacare.com/homecare.

Invacare® Formula™ CG Tilt/Recline for Storm Series® Power Wheelchair

ELECTRONICS Continued

ASL OPTIONAL ACCESSORIES

- EGSBLK Egg Switch-Black \$156.00
- ASL611 Adjustable Switch Mounting Bracket \$252.00

MISCELLANEOUS

BATTERY BOXES

- GP24 Group 24 Battery Boxes No Charge

BATTERIES

- 24BATTERY 24 Gel Battery (On Chair) \$410.00 E2363

BATTERY CHARGER OPTION

- 110CHARGER 110 Volt Battery Charger No Charge
- L8900 Omit Battery Charger Less \$80.00

MISCELLANEOUS ACCESSORIES

- 1303 Wheel Locks \$45.00
- 9205-3 Suspension Shocks \$400.00 E1016
- FDR Fender \$320.00
Not available with TRBKTS.

FRONT RIGGINGS

CENTER PIVOT STYLE (ASB) (PAIR)

- 70NHD 70° Swingaway Footrests 13 1/4" - 17" Standard
- 70STAPER 70° Tapered Footrests 13" - 17" No Charge
Standard with impact guards on the footrests.
- AT5044 Elevating Legrests \$255.00 E0990
Manual elevating legrest.

FOOTPLATE OPTIONS (PAIR)

- 1651 Flip-Up Composite Footplates No Charge
- 1350 Extra Large Aluminum Footplates No Charge
- AT5543 Adjustable Angle Flip Up Footplates \$165.00 K0040

FOOTREST/LEGREST ACCESSORIES

- 1337 Calf Strap \$27.00 E0038
- IMPACTG Impact Guards \$51.00
- HLOOP Heel Loops (pair) \$45.00 E0951
- LPT2 Longer Pivot Slide Tube For 1350 \$83.00
- ALPT Longer Pivot And Slide Tube \$84.00
- ALPT4 Longer Pivot And Slide Tube 4" \$84.00

Conventional Seat Sizes	
Seat Widths	16" - 22"
Seat Depths	16" - 22"
Back Heights	20" - 26"

NOTE: All specifications and dimensions are approximate.
For the most current pricing information, go to www.invacare.com/homecare.

Invacare® Formula™ CG Tilt/Recline for Storm Series® Power Wheelchair

Memory Card - BASIC

- Standard with ALL Expandable Control Systems (CMPJ+, MK6i Display).
- Used to back up / restore programmed settings for only one chair.
- May be given to user or kept with provider's files for safe keeping.
- Does not contain advanced diagnostics / help screens / file structure.
- Not compatible with SPJ+ Joysticks.

Memory Card - PROFESSIONAL

- Standard with ALL MK6i Programmers.
- Available also with USB Card Reader on Order Form.
- Contains advanced diagnostics, help screens, software upgrades, file storage / retrieval
- Used to back up multiple chairs and programming settings.
- Intended for use by qualified providers only.
- Not compatible with SPJ+ Joysticks.

After Market Backrests

The following companies also have hardware to make their backrests work with our new canes:

- AEL - planar seating system
- Freedom Designs - planar seating system
- Canyon - planar seating system

Seating and Positioning

For your seating needs, please refer to the following rehab seating and positioning products: Price List Form No. 07-027

Batteries (not included)

GP24 batteries with terminal configuration (negative on the left and positive on the right) **MUST** be used. GP24 batteries that have the reverse terminal configuration **MUST** not be used. Terminals **MUST** have a cross hole for proper battery connection. See Owner's Manual, part number 1114809. These recommendations **MUST** be followed otherwise injury and damage may occur.

Tech Tips

Operating Powered Seating through the Driver's Control. All chairs require adding or replacing one component to allow actuator operation through the driver control.

Multiple Actuator Systems:

Replace the standard 4 Way Switch box with a S4WSB (Multiple Actuator Interface Box) - E2311.

Formula CG Tilt/Recline Features and Specifications

- Weight capacity of 300 lbs.
- Recline range of 170°.
- Mechanical Shear reduction of 3".
- Conventional or Contoura style seating.
- Adjustable seat frame with Conventional style.
- Seat-To-Floor height of 17 3/4" at 0° tilt.

Comments

NOTE: All specifications and dimensions are approximate.
For the most current pricing information, go to www.invacare.com/homecare.

Invacare Corporation

www.invacare.com

USA

One Invacare Way

Elyria, OH USA

44035-4190

800-333-6900

©2016 Invacare Corporation

All rights reserved. Trademarks are identified by the symbols™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications and prices are subject to change without notification.



Yes, you can.®

Revised December 1, 2016